Ball and Vase by Mr. Magic

The performer shows a **Ball Vase**. The **Ball Vase** is opened, and a ball is seen inside the vase. The ball is withdrawn and placed into the pocket. The empty vase is closed. After a while the vase is reopened. The ball is seen returned from the pocket to the vase! The lid is closed. One hand is then positioned under the vase. The ball penetrates through the vase and drops into the hand! The lid is opened to show that the vase is empty.

As you play around with the apparatus, you will find several routines possible, including all the standard routines that can be performed with a normal Ball and Vase set.

Fun to perform, and very mystifying for your audience!